

PQ001 MRI Contrast Agent Analyzer

PRODUCT DESCRIPTION


The PQ001 NMR Analyzer was launched in 2008. After years of upgrading, PQ001 has many advantages such as small size, high precision, good repeatability, good stability and excellent cost/benefit performance. Based on these advantages, PQ001 has been widely used in MRI Contrast Agent Researches.

Basic parameters:

1. Magnet: permanent magnet;
2. Magnetic field intensity: $0.5 \pm 0.08T$
3. Probe: $\varnothing 15mm$
4. Size (LxWxH): 1685mm×520mm×386mm;
5. Weight: 134Kg;

Functions:

- Relaxation analysis of T2*, T2 and T1;



T1 Relaxation Time

T2 Relaxation Time

