

PQ001 SFC NMR Analyzer Solid Fat Content Analyzer

PRODUCT DESCRIPTION

The PQ001 SFC NMR Analyzer was launched in 2008. After years of upgrading, PQ001 has many advantages such as small size, high precision, good repeatability, good stability and excellent cost/benefit performance. Based on these advantages, PQ001 has been widely used in the determination of Solid Fat Content.

Basic parameters:

1. Magnet: permanent magnet;
2. Magnetic field intensity: $0.5 \pm 0.08T$
3. Probe: $\varnothing 10mm$
4. Size (LxWxH): 1685mm×520mm×386mm;
5. Weight: 134Kg;

Functions:

- Determination of solid fat content(Cocoa butter, Margarine and butter etc.);

