

## PQ001 Spin Finish NMR Analyzer

### PRODUCT DESCRIPTION

PQ001 NMR Analyzer was launched in 2008. After years of upgrading, PQ001 has many advantages such as small size, high precision, good repeatability, good stability and excellent cost/benefit characteristics. Based on these advantages, PQ001 has been widely used in the determination of oil content of textile fiber.

---

### Basic parameters:

1. Magnet: permanent magnet;
2. Magnetic field intensity:  $0.5 \pm 0.08T$
3. Probe:  $\varnothing 25mm$
4. Size (LxWxH): 1685mm×520mm×386mm;
5. Weight: 134Kg;

---

### Functions:

- – Rapid determination of Spin Finish content of textile fiber;;
- 

